1. Indication
GR-14.1 denture is a Methacrylate based resin for 3D-printing systems (SLA, DLP) with irradiation sources (≤ 405 nm) for the production of denture bases.

2. Contraindication
GR-14.1 denture restorations are contraindicated ...

- ... if a patient is known to be allergic to any of the ingredients.
- ... direct intraoral contact of uncured or partially uncured material.
- ... for every application that is not part of the indication (see above).

3. Requirements

1. File Format
 STL (see geometric presettings for CAD)

2. Hardware (3D Printing) – Information available from:
 ASIGA - 2/19-21 Bourke Road, Alexandria, NSW 2015, Australia
 Structo Pte Ltd - 60 Mac Pherson Road, #06-08A, The Siemens Centre, Singapore (348615)
 MIICRAFT - 27 Queen St. East Suite 1401 Toronto, Ontario MSC 2M6, Canada
 DAZZ3D - Rm K-2F, Building A, No. 9 Zhongshanyuan Rd., Shenzhen, China 518052

3. Hardware (Post Curing) – Information available from:
 pro3dure medical (CD-1, CD-2) - Am Burgberg 13, D-58642 Iserlohn, Germany
 Structo Pte Ltd - 60 Mac Pherson Road, #06-08A, The Siemens Centre, Singapore (348615)
 NK Optik (Otoflash G171) - Isarstraße 2, 82065 Baierbrunn, Germany
 Prusa Research a.s. (Prusa CW1) - Partyzánská 188/7a, 170 00 Praha 7, Czech Republic

4. Bonding agent for fixing denture teeth in the denture base – Information available from:
 VITA Zahnfabrik H. Rauter GmbH & Co. KG - Spitalgasse 3, D-79713 Bad Säckingen, Germany

4. Material
GR-14.1 denture consists of functional methacrylic resins, initiators, dyes and stabilisers.

5. Geometric presetting
Lingual Ridge ≥ 2.5 mm (Maxilla / Mandibula)
Palatal/Lingual ≥ 3.0 mm (Maxilla) / ≥ 2.5 (Mandibula)
Facial/Buccal ≥ 2.5 mm (Maxilla / Mandibula)

6. Material parameters
 Radiation penetration depth controllable by exposure time
 ≥ 50 μ m (4 mils)
 ≥ 100 μ m (5 mils)

7. Manufacturing process (fig. 1–10, page 1)


- Prepare data (CAD & build preparation).
- Choose process parameters (build style, etc.).
- Transfer prepared data to 3D printing device.
- Prepare 3D printing – shake bottle.
- Fill resin tank of 3D printing device.
- Build the parts.
- Clean parts (with IPA $\geq 97\%$ or equivalent cleaning agent) approx. 4 min. in an ultrasonic bath or equivalent device – precleaning recommended.
- Dry parts (with compressed air until no residues of IPA or equivalent cleaning agent).
- Post curing (10 min.): inert atmosphere recommended (use adequate light curing device).
- Finish parts.

8. Bonding of denture teeth
 If using printed denture teeth – see additional IFU **GR-17.1 temporary It** for details. In order to realize the optimal bonding of artificial teeth to the dental object, it can be necessary to roughen the surface of the denture and to use an appropriate adhesive. The suitability of a bonding agent or adhesive must be approved prior use. Recommended bonding agent is VITA VIONIC® BOND. The commonly used techniques and dental work steps must be carefully considered when manufacturing the dentures. Please follow the instructions and indications of the corresponding suppliers.

9. Notice
 Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. Ensure personal protective gear during processing.
Caution: Polymerized resins are chemically resistant – avoid stains on clothing! Avoid any contact with skin and eyes! In case of accidental contact, rinse with adequate running water, consult a doctor if necessary. The lot number and the best before date are indicated on each packaging. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date. Dispose of contents/container according to official regulations.

10. Hazard statements
 May cause an allergic skin reaction. May cause long lasting harmful effects to aquatic life.

Instructions for use · Gebrauchsanweisung
 Mode d'emploi · Instrucciones de uso
 Istruzioni per l'uso

Physical properties*/
 Physikalische Eigenschaften*/
 Propriétés physiques*/
 Propiedades físicas*/
 Proprietà fisiche*:

GR-14.1 denture
 ■ Colours/ Farben/ couleurs/
 colores/ colori:

orange-pink/ orange-pink/
 orange-rose/ rosa-naranja/
 rosa-arancione

light-pink/ hell-pink/
 rose-clair/ rosa-claro/
 rosa-chiaro

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1 kg:
 light-pink/ hell-pink/
 rose-clair/ rosa-claro/
 rosa-chiaro
 REF: D1001303

1 kg:
 deep-pink/ dunkel-pink/
 rose-profond/ rosa-
 profundo/ rosa-scuro
 REF: D1001304

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Distributor (US):
 pro3dure medical LLC
 9825 Valley View Road · Eden Prairie, MN 55344
 Phone: 952-426-1928 · Fax: 952-681-7515
 info@pro3dure.com · www.pro3dure.com · Made in Germany

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